

## Instructions for use:

1. Remove all excess blood from the site of the bleeding
2. Immediately apply a liberal amount of HaemoCer™ directly to the source of the bleeding
3. Apply appropriate pressure over the wound by using a dry gauze
4. Before removing the gauze, irrigation is recommended to dissolve adhesions to the gauze
5. Remove excess HaemoCer™ powder by irrigation
6. If bleeding continues, repeat the application of HaemoCer™



Reference Number	Product	Content
HFP105	HaemoCer™ 5g incl. applicator	5 pcs




### BioCER Entwicklungs-GmbH

Ludwig-Thoma-Straße 36c · 95447 Bayreuth

Tel. +49 (0921) 78 77 70-0 · Fax +49 (0921) 78 77 70 79

info@biocer-gmbh.de · [www.biocer-gmbh.de](http://www.biocer-gmbh.de)



HDD10801/2010-12

PLANT BASED  
**100%**  
ABSORBABLE

# BioCER

Effective and safe haemostasis  
derived from nature.  
**MADE IN GERMANY**

# HaemoCer™

## HAEMOSTATIC POWDER

# HaemoCer™ stops bleeding

# FAST – EFFECTIVE - SAFE

## Haemostasis derived from nature

- HaemoCer™ is an adjunct haemostatic agent for use in surgical procedures.
- HaemoCer™ consists to 100% of purified plant based starch and contains no animal or human components.
- HaemoCer™ is applied directly to the source of the bleeding and should cover the wound site completely.
- The hygroscopic particles of HaemoCer™ absorb liquid from the blood and initiate a dehydration process.
- Solid blood components are concentrated and building a clot. The physiological clotting cascade is accelerated without any chemical or pharmaceutical interaction.
- After haemostasis is achieved, a transparent and adhesive gel film protects the wound which also acts as a barrier against further blood discharge.
- HaemoCer™ is completely absorbed within 48h without leaving any residues behind.



## Versatile use

- Treatment of diffuse and local bleedings
- For bleedings around suture lines
- Substitution of fibrin glue for haemostasis in vascular surgery
- Substitution of bone wax and electrocautery cauter for haemostasis of bone bleedings

## Fields of application

- General Surgery
- Cardiac and Vascular Surgery
- Tumour Surgery
- Oral and maxillofacial Surgery
- Plastic Surgery
- NeuroSurgery
- Urology / Gynaecology
- Orthopaedics
- Thoracic Surgery
- Transplant Surgery

**100% absorbable**



- ▶ Effective and safe haemostasis supports the clinical time management
- ▶ Reduction of homologous transfusions and blood donations
- ▶ Decrease of foreign body reactions and complications due to rapid and residue free absorption
- ▶ Ready for use without time-consuming preparations (defrosting, assembling etc.)

- ▶ Multidisciplinary use for surgical bleedings
- ▶ Easy to apply even in difficult bleeding locations
- ▶ Safe, fast and effective without interaction or damage of the surrounding tissue
- ▶ HaemoCer™ is not affected by the patient's clotting status